

Mouse CD276/B7-H3 & CD3E Bispecific, InVivo Block/Neutralizing Antibody

Summary

Catalog No.	BGT-ANT-37451
Species reactivity	Mouse
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Clonality	Monoclonal
Isotype	V-Kappa-VH _ V-Kappa-VH'
Applications	FuncS
Target	CD276, B7-H3, Costimulatory molecule, B7H3, 4Ig-B7-H3, B7 homolog 3, CD276 antigen, T3E, T-cell surface antigen T3/Leu-4 epsilon chain, CD3e, CD3E, T-cell surface glycoprotein CD3 epsilon chain
Purification	Purified by Nickel column.
Endotoxin level	Please contact with the lab for this information.
Accession	Q8VE98 & P22646
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only. Not suitable for clinical or therapeutic use.