


Mouse GIPR/GIP-R, InVivo Block/Neutralizing Antibody

Summary


Catalog No.	BGT-ANT-37455
Host species	Mouse
Species reactivity	Mouse
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Clonality	Monoclonal
Isotype	IgG2a, kappa
Applications	ELISA, Neutralization
Target	Gastric inhibitory polypeptide receptor, GIP-R, Glucose-dependent insulintropic polypeptide receptor, Gipr
Purification	Protein A or G purified from cell culture supernatant.
Endotoxin level	Please contact with the lab for this information.
Accession	Q0P543
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only. Not suitable for clinical or therapeutic use.

Data Image



Detects Mouse GIPR/GIP-R in indirect ELISAs.

Bioactivity



SDS-PAGE