


Mouse IL17A Antibody, InVivo Block/Neutralizing Antibody

Overview


Catalog No.	BGT-ANT-37456
Host species	Mouse
Species reactivity	Mouse
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Purity	>95% as determined by SDS-PAGE.
Clonality	Monoclonal
Isotype	IgG2a
Applications	FCM, Neutralization
Target	CTLA-8, Cytotoxic T-lymphocyte-associated antigen 8, IL17, IL-17A, IL17A, CTLA8, IL-17, Interleukin-17A
Purification	Protein A/G purified from cell culture supernatant.
Endotoxin level	Please contact with the lab for this information.
Accession	Q62386
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Clone ID	BZN035
Note	For research use only. Not suitable for clinical or therapeutic use.

Images




SDS-PAGE for Anti-Mouse IL17A Antibody (BZN035)

SDS-PAGE



Detects Mouse IL17A in indirect ELISAs.

Bioactivity



SEC-HPLC detection for Anti-Mouse IL17A Antibody (BZN035).

Bioactivity