

Pseudorabies Virus/PRV gB/Envelope glycoprotein B/UL27 Antibody, InVivo Block/Neutralizing Antibody

Summary

Catalog No.	BGT-ANT-37468
Clone ID	1H1
Host species	Mouse
Species reactivity	Suid alphaherpesvirus 1
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Clonality	Monoclonal
Isotype	IgG2a, kappa
Applications	ELISA, Neutralization
Target	Envelope glycoprotein B, gB, UL27, gB, UL28
Purification	Protein A/G purified from cell culture supernatant.
Endotoxin level	Please contact with the lab for this information.
Accession	T2FL65
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.