

Human APN Protein (His-Tag) Recombinant Protein

Cat #: BGT-PCT-10141

Size: 50ug 500ug 1mg

Shipping: The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.

Storage: Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.

Synonyms: Aminopeptidase N; ANPEP; AP-M; APN; AP-N; CD13 antigen; CD13; CD13APN; PEPN; PEPNhAPN

Product Information

Source: Human Cells

Description: Recombinant Human Aminopeptidase N is produced by our Mammalian expression system and the target gene encoding Lys69-Lys967 is expressed with a 6His tag at the C-terminus.


Accession: P15144

Predicted Mol Mass: 103.5 KDa

Apparent Mol Mass: 110-130 KDa, reducing conditions

Formulation: Supplied as a 0.2 μm filtered solution of PBS, 5% Trehalose, pH 7.4

Purity-SDS-PAGE: Greater than 95% as determined by reducing SDS-PAGE.



Note:

The product listed herein is for research use only and is not intended for use in human or clinical diagnosis.