

Human PFKM (His-Tag) Recombinant Protein

Cat #: BGT-PCT-11644

Size: 50ug 500ug 1mg

Storage: Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.

Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.

Please minimize freeze-thaw cycles.

Shipping: The product is shipped on dry ice/polar packs.

Upon receipt, store it immediately at the temperature listed below.

Synonyms: 6-phosphofructokinase, muscle type; Phosphofructo-1-kinase isozyme A; Phosphofructokinase 1; Phosphohexokinase; PFKM; PFKX

Product Information

Source: Human Cells

Description: Recombinant Human PhosphoFructoKinase, Muscle Type is produced by our Mammalian expression system and the target gene encoding Thr2-Val780 is expressed with a 6His tag at the C-terminus.

Accession: P08237


Predicted Mol Mass: 86.1 KDa

Apparent Mol Mass: 93 KDa, reducing conditions

Endotoxin: $< 1 \text{ EU}/\mu\text{g}$ as determined by LAL test

Formulation: Supplied as a $0.2 \mu\text{m}$ filtered solution of 20mM PB, 150mM NaCl, 5mM EDTA, 5% Trehalose, pH 6.9.

Purity: Greater than 95% as determined by reducing SDS-PAGE

**Note:**

The product listed herein is for research use only and is not intended for use in human or clinical diagnosis.