

Human FSH (Flag/His-Tag) Recombinant Protein

Cat #: BGT-PCT-11901

Size: 50ug 500ug 1mg

Shipping: Ambient temperature. Upon receipt, store it immediately at the temperature listed below.

Storage: Lyophilized protein should be stored at $\leq -20^{\circ}\text{C}$, stable for one year after receipt. Reconstituted protein solution can be stored at $2-8^{\circ}\text{C}$ for 2-7 days. Aliquots of reconstituted samples are stable at $\leq -20^{\circ}\text{C}$ for 3 months.

Product Information

Source: Human Cells

Description: Recombinant Human Follicle-Stimulating Hormone is produced by our Mammalian expression system and the target gene encoding Ala25-Ser116&Asn19-Glu129 is expressed with a Flag tag&6His at the C-terminus.

Accession: P01215 & P01225

Predicted Mol Mass: 11.5&13.5 KDa


Apparent Mol Mass: 20-30 KDa, reducing conditions

Formulation: Lyophilized from a $0.2\text{ }\mu\text{m}$ filtered solution of PBS, pH 7.4.

Reconstitution: Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than $100\mu\text{g/ml}$. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Reconstitution: Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than $100\mu\text{g/ml}$. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Purity-SDS-PAGE: Greater than 95% as determined by reducing SDS-PAGE.



Note:

The product listed herein is for research use only and is not intended for use in human or clinical diagnosis.