

Biotinylated Human Fibronectin (ED-B domain Avi/His-Tag)

Recombinant Protein

Cat #: BGT-PCT-12093

Size: 20ug 100ug

Storage: Lyophilized protein should be stored at $\leq -20^{\circ}\text{C}$, stable for one year after receipt.

Reconstituted protein solution can be stored at $2-8^{\circ}\text{C}$ for 2-7 days.

Aliquots of reconstituted samples are stable at $\leq -20^{\circ}\text{C}$ for 3 months.

Shipping: The product is shipped at ambient temperature.

Upon receipt, store it immediately at the storage temperature.

Synonyms: Fibronectin; FN1

Product Information

Source: E.coli

Description: Biotinylated Recombinant Human Fibronectin is produced by our E.coli expression system and the target gene encoding Glu1266-Thr1356 is expressed with a Avi, 6His tag at the N-terminus

Accession: P02751-15

Predicted Mol Mass: 13.4 KDa


Apparent Mol Mass: 15 KDa, reducing conditions

Endotoxin: $< 1 \text{ EU}/\mu\text{g}$ as determined by LAL test

Formulation: Lyophilized from a $0.2 \mu\text{m}$ filtered solution of PBS, pH 7.4

Purity: Greater than 95% as determined by reducing SDS-PAGE

Reconstitution: Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than $100\mu\text{g}/\text{ml}$. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

**Note:**

The product listed herein is for research use only and is not intended for use in human or clinical diagnosis.