

Human Granzyme B (C-6His) Recombinant Protein

Cat #: BGT-PCT-12622

Size: 50ug 500ug 1mg

Shipping: The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.

Storage: Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.

Product Information

Source: Human Cells

Description: Recombinant Human Granzyme B is produced by our Mammalian expression system and the target gene encoding Ile21-Tyr247 is expressed with a 6His tag at C terminus. This enzyme has been activated during the purification process.

Accession: P10144


Predicted Mol Mass: 26.6 KDa

Apparent Mol Mass: 32-42 KDa, reducing conditions

Formulation: Supplied as a 0.2 μm filtered solution of PBS, pH7.4.

Endotoxin: $< 1 \text{ EU}/\mu\text{g}$ as determined by LAL test.

Purity-SDS-PAGE: Greater than 95% as determined by reducing SDS-PAGE.



Note:

The product listed herein is for research use only and is not intended for use in human or clinical diagnosis.