

Dog NGF/Beta-NGF Recombinant Protein


Summary

Catalog No.	BGT-PPT-10157
Alternative Names	Beta-nerve growth factor, NGFB, Beta-NGF, NGF
Form	Lyophilized
Storage buffer	Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Purity	>80% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Endotoxin level	Please contact with the lab for this information.
Expression system	Mammalian Cells
Accession	NP_001181879.1
Protein length	Met1-Gly238, C-Fc
Nature	Recombinant
Predicted molecular weight	51.9 kDa
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Reconstitution	Reconstitute in sterile water for a stock solution.
Species	Canis lupus familiaris (Dog) (Canis familiaris)
Shipping	In general, proteins are provided as lyophilized powder/frozen liquid. They are shipped out with dry ice/blue ice unless customers require otherwise.

Note

For research use only.

Data Image



SDS PAGE for recombinant Dog NGF / Beta-NGF

SDS-PAGE