

Human ATP2C1 Recombinant Protein

Summary

Catalog No.	BGT-PPT-10738
Alternative Names	Calcium-transporting ATPase type 2C member 1, ATP-dependent Ca(2+) pump PMR1, PMR1L, KIAA1347, ATP2C1, ATPase 2C1
Form	Lyophilized
Storage buffer	Lyophilized from a solution in PBS pH 7.4, 0.02% N-Lauroylsarcosine, 1 mM EDTA, 4% Trehalose, 1% Mannitol.
Purity	>90% as determined by SDS-PAGE.
Applications	SDS-PAGE, WB, ELISA, Immunogen, Bioactivity testing in progress
Endotoxin level	Please contact with the lab for this information.
Expression system	E. coli
Accession	P98194
Protein length	Val313-Ile699, C-His
Nature	Recombinant
Predicted molecular weight	42.82 kDa
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles.Store at 2 to 8 °C for one week .Store at -20 °C for twelve months from the date of receipt.
Reconstitution	Reconstitute in sterile water for a stock solution.
Species	Human

Shipping

In general, proteins are provided as lyophilized powder/frozen liquid. They are shipped out with dry ice/blue ice unless customers require otherwise.

Note

For research use only.

Data Image



SDS-PAGE for Recombinant Human ATP2C1 protein

SDS-PAGE