

Human CXCL13/BCA-1/BLC Recombinant Protein

Overview

Catalog No.	BGT-PPT-12142
Form	Lyophilized
Storage buffer	Lyophilized from a solution in PBS pH 7.4, 0.02% N-Lauroylsarcosine, 1 mM EDTA, 4% Trehalose, 1% Mannitol.
Purity	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Endotoxin level	Please contact with the lab for this information.
Expression system	E. coli
Accession	O43927
Protein length	Val23-Arg94, N-His
Nature	Recombinant
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Reconstitution	Reconstitute in sterile water for a stock solution. A copy of datasheet will be provided with the products, please refer to it for details.
Alternative name	SCYB13, BLC, CXCL13, B cell-attracting chemokine 1, B lymphocyte chemoattractant, Small-inducible cytokine B13, CXC chemokine BLC, Angie, C-X-C motif chemokine 13, BCA-1, BCA1
Species	Homo sapiens (Human)

Shipping

In general, proteins are provided as lyophilized powder/frozen liquid. They are shipped out with dry ice/blue ice unless customers require otherwise.

Predicted molecular weight

10.98 kDa

Note

For research use only.