

# Human FOXO3 Recombinant Protein

## Summary

---

|                                   |                                                                                                                                                                           |
|-----------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Catalog No.</b>                | BGT-PPT-12768                                                                                                                                                             |
| <b>Alternative Names</b>          | Forkhead in rhabdomyosarcoma-like 1, AF6q21 protein, Forkhead box protein O3, FOXO3, FKHL1, FOXO3A                                                                        |
| <b>Form</b>                       | Lyophilized                                                                                                                                                               |
| <b>Storage buffer</b>             | Lyophilized from a solution in PBS pH 7.4, 0.02% N-Lauroylsarcosine, 1mM EDTA, 4% Trehalose, 1% Mannitol.                                                                 |
| <b>Purity</b>                     | >90% as determined by SDS-PAGE.                                                                                                                                           |
| <b>Applications</b>               | ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress                                                                                                           |
| <b>Endotoxin level</b>            | Please contact with the lab for this information.                                                                                                                         |
| <b>Expression system</b>          | E. coli                                                                                                                                                                   |
| <b>Accession</b>                  | O43524                                                                                                                                                                    |
| <b>Protein length</b>             | Met372-Gly673, N-His                                                                                                                                                      |
| <b>Nature</b>                     | Recombinant                                                                                                                                                               |
| <b>Predicted molecular weight</b> | 34.46 kDa                                                                                                                                                                 |
| <b>Stability and Storage</b>      | Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt. |
| <b>Reconstitution</b>             | Reconstitute in sterile water for a stock solution. A copy of datasheet will be provided with the products, please refer to it for details.                               |
| <b>Species</b>                    | Homo sapiens (Human)                                                                                                                                                      |

**Shipping**


In general, proteins are provided as lyophilized powder/frozen liquid. They are shipped out with dry ice/blue ice unless customers require otherwise.

**Note**

For research use only.

## Data Image

---



SDS-PAGE

SDS PAGE for recombinant Human FOXO3 protein