

Human LAMA1 Recombinant Protein

Summary


Catalog No.	BGT-PPT-13602
Alternative Names	Laminin subunit alpha-1, LAMA, S-LAM alpha, S-laminin subunit alpha, Laminin-3 subunit alpha, LAMA1, Laminin A chain, Laminin-1 subunit alpha, laminin 1, laminin 111, EHS laminin, laminin 3, laminin 121, S laminin
Form	Lyophilized
Storage buffer	Lyophilized from a solution in PBS pH 7.4, 0.02% N-Lauroylsarcosine, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Purity	>90% as determined by SDS-PAGE.
Applications	SDS-PAGE,WB,ELISA,Immunogen,Bioactivity testing in progress
Endotoxin level	Please contact with the lab for this information.
Expression system	E. coli
Accession	P25391
Protein length	Val2613-Phe2746, N-His
Nature	Recombinant
Predicted molecular weight	17.25 kDa
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8 °C for one week. Store at -20 to -80 °C for twelve months from the date of receipt.
Reconstitution	Reconstitute in sterile water for a stock solution. A copy of datasheet will be provided with the products, please refer to it for details.

Species Human

Shipping In general, proteins are provided as lyophilized powder/frozen liquid. They are shipped out with dry ice/blue ice unless customers require otherwise.

Note For research use only.

Data Image



SDS-PAGE

SDS PAGE for recombinant Human LAMA1