

Human MUC5AC Recombinant Protein

Overview

Catalog No.	BGT-PPT-14048
Form	Lyophilized
Storage buffer	Lyophilized from a solution in PBS pH 7.4, 1 mM EDTA, 4% Trehalose, 1% Mannitol.
Purity	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Endotoxin level	Please contact with the lab for this information.
Expression system	E. coli
Accession	P98088
Protein length	Tyr5568-His5654, N-His-SUMO
Nature	Recombinant
Predicted molecular weight	21.97 kDa
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Reconstitution	Reconstitute in sterile water for a stock solution. A copy of datasheet will be provided with the products, please refer to it for details.
Alternative name	Mucin-5 subtype AC, tracheobronchial, Gastric mucin, Mucin-5AC, TBM, MUC-5AC, Tracheobronchial mucin, Major airway glycoprotein, MUC5AC, MUC5
Species	Homo sapiens (Human)


Shipping

In general, proteins are provided as lyophilized powder/frozen liquid. They are shipped out with dry ice/blue ice unless customers require otherwise.

Note

For research use only.

Images



SDS-PAGE for Recombinant Human MUC5AC Protein, N-His-SUMO

SDS-PAGE