

Human NCF1 Recombinant Protein

Summary

Catalog No.	BGT-PPT-14126
Alternative Names	SH3PXD1A, NCF1, p47-phox, 47 kDa autosomal chronic granulomatous disease protein, Nox-organizing protein 2, NOXO2, SH3 and PX domain-containing protein 1A, NCF-1, 47 kDa neutrophil oxidase factor, NCF-47K, Nox organizer 2, Neutrophil NADPH oxidase factor 1, Neutrophil cytosol factor 1
Form	Lyophilized
Storage buffer	Lyophilized from a solution in PBS pH 7.4, 0.02% N-Lauroylsarcosine, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Purity	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Endotoxin level	Please contact with the lab for this information.
Expression system	E. coli
Accession	P14598
Protein length	Met1-Val390, N-His
Nature	Recombinant
Predicted molecular weight	46.99 kDa
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Reconstitution	Reconstitute in sterile water for a stock solution. A copy of datasheet will be provided with the products, please refer to it for details.

Species Homo sapiens (Human)


Shipping

In general, proteins are provided as lyophilized powder/frozen liquid.
They are shipped out with dry ice/blue ice unless customers require otherwise.

Note

For research use only.

Data Image



SDS PAGE for recombinant Human NCF1

SDS-PAGE