

Human NFE2L2 Recombinant Protein

Summary

Catalog No.	BGT-PPT-14185
Alternative Names	NFE2-related factor 2, Nuclear factor, erythroid derived 2, like 2, NRF2, Nuclear factor erythroid 2-related factor 2, NFE2L2, Nrf-2, HEBP1, NF-E2-related factor 2
Form	Lyophilized
Storage buffer	Lyophilized from a solution in PBS pH 7.4, 0.02%N-Lauroylsarcosine, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Purity	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Endotoxin level	Please contact with the lab for this information.
Expression system	E. coli
Accession	Q16236
Protein length	Met17-Asn605, N-His
Nature	Recombinant
Predicted molecular weight	68.39 kDa
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Reconstitution	Reconstitute in sterile water for a stock solution. A copy of datasheet will be provided with the products, please refer to it for details.
Species	Homo sapiens (Human)

Shipping

In general, proteins are provided as lyophilized powder/frozen liquid. They are shipped out with dry ice/blue ice unless customers require otherwise.

Note

For research use only.

Data Image

