

# Human NPC1L1 Recombinant Protein

## Summary

---

|                                   |                                                                                                                                                                        |
|-----------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Catalog No.</b>                | BGT-PPT-14253                                                                                                                                                          |
| <b>Alternative Names</b>          | NPC1L1, Niemann-Pick C1-like protein 1, NPC1-like intracellular cholesterol transporter 1                                                                              |
| <b>Form</b>                       | Lyophilized                                                                                                                                                            |
| <b>Storage buffer</b>             | Lyophilized from a solution in PBS pH 7.4, 0.02% N-Lauroylsarcosine, 1mM EDTA, 4% Trehalose, 1% Mannitol.                                                              |
| <b>Purity</b>                     | >90% as determined by SDS-PAGE.                                                                                                                                        |
| <b>Applications</b>               | SDS-PAGE, WB, ELISA, Immunogen, Bioactivity testing in progress                                                                                                        |
| <b>Endotoxin level</b>            | Please contact with the lab for this information.                                                                                                                      |
| <b>Expression system</b>          | E. coli                                                                                                                                                                |
| <b>Accession</b>                  | Q9UHC9                                                                                                                                                                 |
| <b>Protein length</b>             | Arg450-Ser620, N-His                                                                                                                                                   |
| <b>Nature</b>                     | Recombinant                                                                                                                                                            |
| <b>Predicted molecular weight</b> | 21.74 kDa                                                                                                                                                              |
| <b>Stability and Storage</b>      | Use a manual defrost freezer and avoid repeated freeze thaw cycles.Store at 2 to 8 °C for one week .Store at -20 to -80 °C for twelve months from the date of receipt. |
| <b>Reconstitution</b>             | Reconstitute in sterile water for a stock solution.A copy of datasheet will be provided with the products, please refer to it for details.                             |
| <b>Species</b>                    | Human                                                                                                                                                                  |

**Shipping**


In general, proteins are provided as lyophilized powder/frozen liquid. They are shipped out with dry ice/blue ice unless customers require otherwise.

**Note**

For research use only.

## Data Image

---



SDS PAGE for Recombinant Human NPC1L1 Protein

SDS-PAGE