


Human PRKRA Recombinant Protein

Summary

Catalog No.	BGT-PPT-14752
Alternative Names	PKR-associated protein X, Protein activator of the interferon-induced protein kinase, PACT, PRKRA, RAX, PKR-associating protein X, Interferon-inducible double-stranded RNA-dependent protein kinase activator A, Protein kinase, interferon-inducible double-stranded RNA-dependent activator
Form	Lyophilized
Storage buffer	Lyophilized from a solution in PBS pH 7.4, 0.02% N-Lauroylsarcosine, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Purity	>90% as determined by SDS-PAGE.
Applications	SDS-PAGE, WB, ELISA, Immunogen, Bioactivity testing in progress
Endotoxin level	Please contact with the lab for this information.
Expression system	E. coli
Accession	O75569
Protein length	Asn126-Lys313, N-His
Nature	Recombinant
Predicted molecular weight	23.00 kDa
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8 °C for one week. Store at -20 to -80 °C for twelve months from the date of receipt.
Reconstitution	Reconstitute in sterile water for a stock solution. A copy of datasheet will be provided with the products, please refer to it for details.

Species	Human
Shipping	In general, proteins are provided as lyophilized powder/frozen liquid. They are shipped out with dry ice/blue ice unless customers require otherwise.
Note	For research use only.

Data Image



SDS PAGE for Recombinant Human PRKRA Protein

SDS-PAGE