

Human PRMT5 Recombinant Protein

Overview


Catalog No.	BGT-PPT-14758
Form	Lyophilized
Storage buffer	Lyophilized from a solution in PBS pH 7.4, 0.02% N-Lauroylsarcosine, 1 mM EDTA, 4% Trehalose, 1% Mannitol.
Purity	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Endotoxin level	Please contact with the lab for this information.
Expression system	E. coli
Accession	O14744
Protein length	Asn295-Gly466, N-His
Nature	Recombinant
Predicted molecular weight	21.71 kDa
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Reconstitution	Reconstitute in sterile water for a stock solution. A copy of datasheet will be provided with the products, please refer to it for details.
Alternative name	PRMT5, SKB1 homolog, 72 kDa ICLn-binding protein, Jak-binding protein 1, Histone-arginine N-methyltransferase PRMT5, SKB1Hs, JBP1, Shk1 kinase-binding protein 1 homolog, SKB1, HRMT1L5, IBP72, Protein arginine N-methyltransferase 5

Species Homo sapiens (Human)

Shipping In general, proteins are provided as lyophilized powder/frozen liquid. They are shipped out with dry ice/blue ice unless customers require otherwise.

Note For research use only.

Images



SDS-PAGE for Recombinant Human PRMT5 Protein,
N-His

SDS-PAGE