


Human RORC/NR1F3/RORgT Recombinant Protein

Summary

Catalog No.	BGT-PPT-15046
Alternative Names	RAR-related orphan receptor C, RORC, NR1F3, Nuclear receptor ROR-gamma, Nuclear receptor RZR-gamma, Retinoid-related orphan receptor-gamma, RZRG, Nuclear receptor subfamily 1 group F member 3, RORG, RORgT
Form	Lyophilized
Storage buffer	Lyophilized from a solution in PBS pH 7.4, 0.02% N-Lauroylsarcosine, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Purity	>85% as determined by SDS-PAGE.
Applications	SDS-PAGE,WB,ELISA,Immunogen,Bioactivity testing in progress
Endotoxin level	Please contact with the lab for this information.
Expression system	E. coli
Accession	B6ZGS6
Protein length	Thr268-Phe494, N-His-SUMO
Nature	Recombinant
Predicted molecular weight	39.11 kDa
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8 °C for one week. Store at -20 to -80 °C for twelve months from the date of receipt.

Reconstitution	Reconstitute in sterile water for a stock solution. A copy of datasheet will be provided with the products, please refer to it for details.
Species	Human
Shipping	In general, proteins are provided as lyophilized powder/frozen liquid. They are shipped out with dry ice/blue ice unless customers require otherwise.
Note	For research use only.

Data Image



SDS PAGE for recombinant Human RORC / NR1F3

SDS-PAGE