


Human RUNX2 Recombinant Protein

Summary

Catalog No.	BGT-PPT-15091
Alternative Names	AML3, OSF2, RUNX2, Core-binding factor subunit alpha-1, PEA2-alpha A, CBF-alpha-1, Osteoblast-specific transcription factor 2, Oncogene AML-3, CBFA1, Runt-related transcription factor 2, PEBP2A, OSF-2, SL3-3 enhancer factor 1 alpha A subunit, PEBP2-alpha A, Acute myeloid leukemia 3 protein, Polyomavirus enhancer-binding protein 2 alpha A subunit, SL3/AKV core-binding factor alpha A subunit
Form	Lyophilized
Storage buffer	Lyophilized from a solution in PBS pH 7.4, 0.02% N-Lauroylsarcosine, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Purity	>90% as determined by SDS-PAGE.
Applications	SDS-PAGE,WB,ELISA,Immunogen,Bioactivity testing in progress
Endotoxin level	Please contact with the lab for this information.
Expression system	E. coli
Accession	Q13950
Protein length	Val219-Val404, N-His
Nature	Recombinant
Predicted molecular weight	22.80 kDa
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles.Store at 2 to 8 °C for one week .Store at -20 to -80 °C for twelve months from the date of receipt.

Reconstitution	Reconstitute in sterile water for a stock solution.A copy of datasheet will be provided with the products, please refer to it for details.
Species	Human
Shipping	In general, proteins are provided as lyophilized powder/frozen liquid. They are shipped out with dry ice/blue ice unless customers require otherwise.
Note	For research use only.

Data Image



SDS PAGE for recombinant Human RUNX2 protein

SDS-PAGE