


Human TIA1 Recombinant Protein

Summary

Catalog No.	BGT-PPT-15670
Alternative Names	TIA1, T-cell-restricted intracellular antigen-1, RNA-binding protein TIA-1, TIA-1, Nucleolysin TIA-1 isoform p40, p40-TIA-1
Form	Lyophilized
Storage buffer	Lyophilized from a solution in PBS pH 7.4, 0.02% N-Lauroylsarcosine, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Purity	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Endotoxin level	Please contact with the lab for this information.
Expression system	E. coli
Accession	P31483
Protein length	Asp104-Trp181, N-His-SUMO & C-Strep
Nature	Recombinant
Predicted molecular weight	22.37 kDa
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Reconstitution	Reconstitute in sterile water for a stock solution. A copy of datasheet will be provided with the products, please refer to it for details.

Species	Homo sapiens (Human)
Shipping	In general, proteins are provided as lyophilized powder/frozen liquid. They are shipped out with dry ice/blue ice unless customers require otherwise.
Note	For research use only.

Data Image



SDS PAGE for Recombinant Human TIA1 Protein

SDS-PAGE