

Human TRIM33/TIF1- γ Recombinant Protein

Summary


Catalog No.	BGT-PPT-15828
Alternative Names	TRIM33, RFG7, KIAA1113, RET-fused gene 7 protein, Tripartite motif-containing protein 33, RING-type E3 ubiquitin transferase TRIM33, E3 ubiquitin-protein ligase TRIM33, Protein Rfg7, TIF1-gamma, Ectodermin homolog, TIF1G, Transcription intermediary factor 1-gamma
Form	Lyophilized
Storage buffer	Lyophilized from a solution in PBS pH 7.4, 0.02% N-Lauroylsarcosine, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Purity	>90% as determined by SDS-PAGE.
Applications	SDS-PAGE, WB, ELISA, Immunogen, Bioactivity testing in progress
Endotoxin level	Please contact with the lab for this information.
Expression system	E. coli
Accession	Q9UPN9
Protein length	Asp882-Ala1087, N-His
Nature	Recombinant
Predicted molecular weight	26.01 kDa
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles.Store at 2 to 8 °C for one week .Store at -20 to -80 °C for twelve months from the date of receipt.
Reconstitution	Reconstitute in sterile water for a stock solution.A copy of datasheet will be provided with the products, please refer to it for details.

Species Human

Shipping In general, proteins are provided as lyophilized powder/frozen liquid. They are shipped out with dry ice/blue ice unless customers require otherwise.

Note For research use only.

Data Image



SDS PAGE for Recombinant Human TRIM33 Protein

SDS-PAGE