

Human TUBB/Tubulin beta Recombinant Protein


Summary

Catalog No.	BGT-PPT-15874
Alternative Names	Tubulin beta-5 chain, TUBB, Tubulin beta chain, TUBB5
Form	Lyophilized
Storage buffer	Lyophilized from a solution in PBS pH 7.4, 0.02% N-Lauroylsarcosine, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Purity	>90% as determined by SDS-PAGE.
Applications	SDS-PAGE, WB, ELISA, Immunogen, Bioactivity testing in progress
Endotoxin level	Please contact with the lab for this information.
Expression system	E. coli
Accession	P07437
Protein length	Val170-Val419, N-His
Nature	Recombinant
Predicted molecular weight	30.67 kDa
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8 °C for one week .Store at -20 to -80 °C for twelve months from the date of receipt.
Reconstitution	Reconstitute in sterile water for a stock solution.
Species	Human
Shipping	In general, proteins are provided as lyophilized powder/frozen liquid. They are shipped out with dry ice/blue ice unless customers require otherwise.

Note

For research use only.

Data Image



SDS PAGE for recombinant Human TUBB protein

SDS-PAGE