

Human VTCN1 Recombinant Protein

Summary

Catalog No.	BGT-PPT-16023
Alternative Names	V-set domain-containing T-cell activation inhibitor 1, B7h.5, B7 homolog 4, B7-H4, T-cell costimulatory molecule B7x, B7H4, VTCN1, Protein B7S1, Immune costimulatory protein B7-H4
Form	Lyophilized
Storage buffer	Lyophilized from a solution in PBS pH 7.4, 0.02% N-Lauroylsarcosine, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Purity	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Endotoxin level	Please contact with the lab for this information.
Expression system	E. coli
Accession	Q7Z7D3
Protein length	Arg34-Thr147, N-His
Nature	Recombinant
Predicted molecular weight	14.82 kDa
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Reconstitution	Reconstitute in sterile water for a stock solution. A copy of datasheet will be provided with the products, please refer to it for details.
Species	Homo sapiens (Human)


Shipping

In general, proteins are provided as lyophilized powder/frozen liquid. They are shipped out with dry ice/blue ice unless customers require otherwise.

Note

For research use only.

Data Image



SDS-PAGE

SDS PAGE for recombinant Human VTCN1 protein