

Human ZBP1 Recombinant Protein

Summary

Catalog No.	BGT-PPT-16103
Alternative Names	DNA-dependent activator of IFN-regulatory factors, DLM1, ZBP1, C20orf183, Z-DNA-binding protein 1, DAI, Tumor stroma and activated macrophage protein DLM-1
Form	Lyophilized
Storage buffer	Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Purity	>90% as determined by SDS-PAGE.
Applications	SDS-PAGE,WB,ELISA,Immunogen,Bioactivity testing in progress
Endotoxin level	Please contact with the lab for this information.
Expression system	E. coli
Accession	Q9H171
Protein length	Arg10-Pro167, N-His
Nature	Recombinant
Predicted molecular weight	20.20 kDa
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles.Store at 2 to 8 °C for one week .Store at -20 to -80 °C for twelve months from the date of receipt.
Reconstitution	Reconstitute in sterile water for a stock solution.
Species	Human


Shipping

In general, proteins are provided as lyophilized powder/frozen liquid. They are shipped out with dry ice/blue ice unless customers require otherwise.

Note

For research use only.

Data Image



SDS PAGE for recombinant Human ZBP1 protein

SDS-PAGE