

Mouse CD254/RANKL/TNFSF11 Recombinant Protein

Summary

Catalog No.	BGT-PPT-16367
Alternative Names	ODF, TNF-related activation-induced cytokine, Osteoclast differentiation factor, TNFSF11, TRANCE, Osteoprotegerin ligand, RANKL, Tumor necrosis factor ligand superfamily member 11, Receptor activator of nuclear factor kappa-B ligand, CD254, OPGL
Form	Lyophilized
Storage buffer	Lyophilized from a solution in PBS pH 7.4, 0.02% N-Lauroylsarcosine, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Purity	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Endotoxin level	Please contact with the lab for this information.
Expression system	E. coli
Accession	O35235
Protein length	Arg72-Asp316, N-His-SUMO
Nature	Recombinant
Predicted molecular weight	29.85 kDa
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Reconstitution

Reconstitute in sterile water for a stock solution. A copy of datasheet will be provided with the products, please refer to it for details.

Species

Mus musculus (Mouse)


Shipping

In general, proteins are provided as lyophilized powder/frozen liquid. They are shipped out with dry ice/blue ice unless customers require otherwise.

Note

For research use only.

Data Image



SDS-PAGE for Recombinant Mouse CD254/RANKL/TNFSF 11 protein

SDS-PAGE