

Mouse TNC/Tenascin Recombinant Protein

Summary

Catalog No.	BGT-PPT-16762
Alternative Names	Myotendinous antigen, Neuronectin, GMEM, Tenascin, Tenascin-C, TN-C, TNC, TN, Cytotactin, GP 150-225, Glioma-associated-extracellular matrix antigen, HXB, Hexabrachion, JI
Form	Lyophilized
Storage buffer	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Purity	>90% as determined by SDS-PAGE.
Applications	SDS-PAGE, WB, ELISA, Immunogen, Bioactivity testing in progress
Endotoxin level	Please contact with the lab for this information.
Expression system	E. coli
Accession	Q80YX1
Protein length	Pro1888-Ala2110, N-His Tag
Nature	Recombinant
Predicted molecular weight	28.01 kDa
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8 °C for one week .Store at -20 °C for twelve months from the date of receipt.
Reconstitution	Reconstitute in sterile water for a stock solution.
Species	Mouse


Shipping

In general, proteins are provided as lyophilized powder/frozen liquid. They are shipped out with dry ice/blue ice unless customers require otherwise.

Note

For research use only.

Data Image



SDS-PAGE

SDS-PAGE for Recombinant Mouse TNC/Tenascin protein