

Native Goat IgG

Summary

Catalog No.	BGT-PPT-16804
Alternative Names	IgG, Immunoglobulin gamma
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Purity	>90% as determined by SDS-PAGE.
Applications	SDS-PAGE, WB, ELISA, Immunogen
Endotoxin level	Please contact with the lab for this information.
Expression system	Native
Protein length	Full length.
Nature	Native
Predicted molecular weight	150 kDa
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8 °C for one week .Store at -20 to -80 °C for twelve months from the date of receipt.
Reconstitution	Reconstitute in sterile water for a stock solution. A copy of datasheet will be provided with the products, please refer to it for details.
Species	Goat
Shipping	In general, proteins are provided as lyophilized powder/frozen liquid. They are shipped out with dry ice/blue ice unless customers require otherwise.
Note	For research use only.