


Stable Receptor Related Protein (RAP) LRP1 Inhibitor

Summary

Catalog No.	BGT-PPT-17009
Alternative Names	Stable Receptor Related Protein (RAP) LRP1 Inhibitor, RAP, LRP1 Inhibitor
Form	Lyophilized
Storage buffer	Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Purity	>90% as determined by SDS-PAGE.
Applications	Inhibitor
Endotoxin level	Please contact with the lab for this information.
Nature	Recombinant
Predicted molecular weight	42.00 kDa
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Reconstitution	Reconstitute in sterile water for a stock solution.
Species	Homo sapiens (Human)
Shipping	In general, proteins are provided as lyophilized powder/frozen liquid. They are shipped out with dry ice/blue ice unless customers require otherwise.
Note	For research use only.

Data Image



SDS-PAGE for Recombinant Stable Receptor
Related Protein (RAP) LRP1 Inhibitor

SDS-PAGE