

EDTA Anticoagulant Cynomolgus Monkey Blood

Cat #: A-CSH029

Size: 10mL

Storage: 2-8 °C (1 month)

Product Description

EDTA Anticoagulant Cynomolgus Monkey Blood is prepared by aseptically collecting cynomolgus monkey blood and treating it with EDTA anticoagulant.

Usage

It should be used according to the specific experiment requirements.

Notes

1. Avoid freeze-thaw cycles.
2. Gently invert and mix daily.
3. Handle aseptically, store at low temperature, and use promptly after opening.
4. Avoid prolonged or long-distance transportation.

Disclaimer

The reagent is only used in the field of scientific research, not suitable for clinical diagnosis or other purposes.