

# Purified Anti-PtdIns(3,4)P2 Antibody

## Overview

---

<b>Catalog No.</b>	BGT-ANT-03175	
<b>Host species</b>	Mouse	
<b>Form</b>	Liquid	
<b>Storage buffer</b>	0.01M PBS pH 7.4	
<b>Purity</b>	>95% as determined by SDS-PAGE.	
<b>Clonality</b>	Monoclonal	
<b>Isotype</b>	IgG2a, kappa	
<b>Applications</b>	ELISA, IB, IF	
<b>Target</b>	PtdIns(3,4)P2; PI-3,4-P2; PIP2	
<b>Purification</b>	Protein A/G purified from cell culture supernatant.	
<b>Endotoxin level</b>	Please contact with the lab for this information.	
<b>Stability and Storage</b>	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.	
<b>Clone ID</b>	SAA2519	
<b>Product Usage Information</b>	Application	Dilution
	ELISA	1:2000-1:20000
<b>Note</b>	For research use only	