

# Human $\alpha$ -MG/Alpha-1-microglobulin Monoclonal Antibody

## Overview

---

<b>Catalog No.</b>	BGT-ANT-03202
<b>Host species</b>	Mouse
<b>Species reactivity</b>	Human
<b>Immunogen</b>	Human $\alpha$ -MG/Alpha-1-microglobulin Protein.
<b>Form</b>	Liquid
<b>Storage buffer</b>	0.01M PBS, pH 7.4.
<b>Clonality</b>	Monoclonal
<b>Isotype</b>	IgG
<b>Applications</b>	ICA, TIA
<b>Target</b>	Protein AMBP, Protein HC, Alpha-1-microglobulin, 1.6.2.-, Alpha-1 microglycoprotein, Complex-forming glycoprotein heterogeneous in charge, Inter-alpha-trypsin inhibitor light chain, ITI-LC, Bikunin, EDC1, HI-30, Uronic-acid-rich protein, Trypstatin, AMBP, HCP, ITIL
<b>Purification</b>	Protein A/G purified from cell culture supernatant.
<b>Stability and Storage</b>	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
<b>Clone ID</b>	1A397
<b>Note</b>	For research use only.