

PEG Recombinant Antibody, clone 6.3

Overview

Catalog No.	BGT-ANT-03205	
Host species	Mouse	
Form	Liquid	
Storage buffer	0.01M PBS, pH 7.4.	
Purity	>95% as determined by SDS-PAGE.	
Clonality	Monoclonal	
Isotype	IgG1, kappa	
Applications	ELISA, FCM, WB	
Target	PEG, Poly(ethylene glycol), Polyethylene glycol, 6.3, 6-3	
Purification	Protein A/G purified from cell culture supernatant.	
Endotoxin level	Please contact with the lab for this information.	
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.	
Clone ID	6.3	
	Application	Dilution
Product Usage Information	ELISA	1:2000-1:20000
	WB	1:500-1:2000
Note	For research use only.	