

Stonin-2 Monoclonal Antibody

Summary

Catalog No.	BGT-ANT-23921
Alternative Names	Stonin-2, STN2, STON2, Stoned B, STNB
Host species	Mouse
Species reactivity	Human, Mouse, Rat
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% by SDS-PAGE.
Clonality	Monoclonal
Isotype	IgG2a, kappa
Applications	IHC, WB
Purification	Protein A or G purified from cell culture supernatant.
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles.Store at +4°C short term (1-2 weeks).Store at -20 °C 12 months. Store at -80°C long term.
Note	For research use only.