

Recombinant Mouse, IL27 EBI3 / Interleukin-27 EBI3,

Animal Free & Carrier Free

Cat #: C-CYP165

Size: 5µg, 20µg, 100µg, 500µg, 1mg

Shipping: Blue Ice

Product Overview

IL-27 protein is a member of the IL-6 superfamily, which is expressed on monocytes, endothelial cells and dendritic cells. It is also referred to as the IL-12 p35-related protein. IL-27 protein is an early product of activated antigen-presenting cells and drives rapid clonal expansion of naive CD4(+) T cells and plays a role in the early regulation of Th1 cells initiation which drives efficient adaptive immune response. It has an antiproliferative activity on melanomas through WSX-1/STAT1 signaling. Thus, IL-27 protein may be an attractive candidate as an antitumor agent applicable to cancer immunotherapy

Product Information

Source: *Escherichia Coli*.

Purity: >95% as determined by SDS-PAGE. Ni-NTA chromatography.

Biological Activity: Measure by its ability to protect HepG2 cells infected with encephalomyocarditis (EMC) virus.


The ED₅₀ for this effect is < 5 ng/mL.

Endotoxin: <0.1 EU per 1 µg of the protein by the LAL method.

Amino acid sequence:

MALVALSQPRVQCHASRYPVAVDCSWTPLQAPNSTRSTSFYATYRLGVATQQSQPCLQRSPQASRCTIPDVHLFSTEBI3VPYMLNVTAV
HPGGASSLLAFVAERIIKPDPPPEGVRLRTAGQRLQVLWHPASWPFDPDIFSLKYRLRYRRRGASHFRQVGPIEATTFTLRNSKPHAKYCIQV
SAQDLTDYGKPSDWSLPGQVESAPHKP with polyhistidine tag at the C-terminus.

Formulation: Lyophilized from a sterile filtered aqueous solution in 1×PBS, pH 4.5.



SDS-PAGE analysis of recombinant mouse IL-27 EB13

Usage Method

1. Before opening, it is recommended to centrifuge at 3000-3500 rpm for 5 minutes.
2. Reconstitute to a concentration of 0.1-1.0 mg/mL in sterile distilled H₂O. Allow the solution to sit at room temperature for at least 20 minutes to ensure complete dissolution. Avoid vigorous vortexing.
3. The reconstituted solution can be stored at 2-8°C for up to 1 week.
4. For long-term storage, it is recommended to further dilute the solution with a carrier protein (such as 0.1% BSA, 10% FBS, or 5% HSA) to a concentration of no less than 10 µg/mL and aliquot for storage at -20°C to -80°C for 3 to 6 months. If serum-free experiments are required, a 5% trehalose solution can be used as a carrier instead. Avoid repeated freeze-thaw cycles.

Storage

Physical Appearance	Storage	Stability
Lyophilized powder	-20°C to -80°C	1 year
Reconstitution (initial)	2°C to 8°C	Less than 1 week
Reconstitution (after dilution)	-20°C to -80°C	3 to 6 months

Note

The reagent is only used in the field of scientific research, not suitable for clinical diagnosis or other purposes.