

Recombinant Human IFN α 1a / Interferon alpha 1a,

Animal Free & Carrier Free

Cat #: C-CYP277

Size: 5 μ g, 20 μ g, 100 μ g, 500 μ g, 1mg

Shipping: Blue Ice

Product Overview

Interferon alpha-1 or interferon alpha -D, encoded by the human IFNA1 gene, is a member of the type I interferon family. It is expressed in peripheral blood leukocytes and lymphoblastoid cells. There are four highly conserved cysteine residues which form two disulfide bonds. One bond is a necessity for biological activity, which plays a critical role in conducting non-specific resistance against a variety of viral infections, and regulation of cell proliferation and modulates immune responses.

Product Information

Source: *Escherichia Coli*.

Purity: >95% as determined by SDS-PAGE. Ni-NTA chromatography.


Biological Activity: Measure by its ability to inhibit IL-8 secretion in human PBMCs in the presence of LPS. The ED₅₀ for this effect is < 1.12 μ g /mL

Endotoxin: <0.1EU per 1 μ g of the protein by the LAL method.

Amino acid sequence:

CDLPETHSLDNRRTLMLLAQMSRISPSSCLMDRHDGFGFPQEEFDGNQFQKAPAISVLHELIQQIFNLFTTKDSSAAWDEDLLDKFCTELYQ
QLNDLEACVMQEERVGETPLMNADSILAVKKYFRRITLYLT EKKYSPCAWEVVRAEIMRSLSLSTNLQERLRRKE with polyhistidine tag at the N-terminus.

Formulation: Lyophilized from a sterile filtered aqueous solution in 1X PBS, pH 8.0



SDS-PAGE analysis of recombinant human IFN alpha 1a

Usage Method

1. Before opening, it is recommended to centrifuge at 3000-3500 rpm for 5 minutes.
2. Reconstitute to a concentration of 0.1-1.0 mg/mL in sterile distilled H₂O. Allow the solution to sit at room temperature for at least 20 minutes to ensure complete dissolution. Avoid vigorous vortexing.
3. The reconstituted solution can be stored at 2-8°C for up to 1 week.
4. For long-term storage, it is recommended to further dilute the solution with a carrier protein (such as 0.1% BSA, 10% FBS, or 5% HSA) to a concentration of no less than 10 µg/mL and aliquot for storage at -20°C to -80°C for 3 to 6 months. If serum-free experiments are required, a 5% trehalose solution can be used as a carrier instead. Avoid repeated freeze-thaw cycles.

Storage

Physical Appearance	Storage	Stability
Lyophilized powder	-20°C to -80°C	1 year
Reconstitution (initial)	2°C to 8°C	Less than 1 week
Reconstitution (after dilution)	-20°C to -80°C	3 to 6 months

Note

The reagent is only used in the field of scientific research, not suitable for clinical diagnosis or other purposes.