

Recombinant Human CNTF / Ciliary Neurotrophic Factor,

Animal Free & Carrier Free

Cat #: C-CYP352

Size: 5µg, 20µg, 100µg, 500µg, 1mg

Shipping: Blue Ice

Product Overview

The ciliary neurotrophic factor is a protein that in humans is encoded by the CNTF gene. It is a hypothalamic neuropeptide that is a potent survival factor for neurons and oligodendrocytes and may be relevant in reducing tissue destruction during inflammatory attacks. CNTF has also been shown to be expressed by cells on the bone surface and to reduce the activity of bone-forming cells (osteoblasts).

Product Information

Source: *Escherichia Coli*.

Purity: >98% as determined by SDS-PAGE. Ni-NTA chromatography.


Biological Activity: Measure by its ability to induce proliferation in TF-1 cells. The ED50 for this effect is <0.15 µg/mL.

Endotoxin: <0.1EU per 1µg of the protein by the LAL method.

Amino acid sequence:

MAFTEHSPLTPHRRDLCSRSIWLARKIRSDLTALTQSYVKHQGLNKNINLDSADGMPVASTDQWSQLTQAQRLQQNLQAYRTFHVLLARL
LQDQQVHFTPTQGDFHQAIHTLLLQVAAFAYQIQQLMILLQYKIPRNQADGMPINVGDDGLFQKKLWGLKVLQQLSQWTVRSIHDLRFI
SSHQTGIPARGSHYIANNKKM with polyhistidine tag at the C-terminus.

Formulation: Lyophilized from a sterile filtered aqueous solution in 1×PBS, pH 7.4.



SDS-PAGE analysis of recombinant human CNTF

Usage Method

1. Before opening, it is recommended to centrifuge at 3000-3500 rpm for 5 minutes.
2. Reconstitute to a concentration of 0.1-1.0 mg/mL in sterile distilled H₂O. Allow the solution to sit at room temperature for at least 20 minutes to ensure complete dissolution. Avoid vigorous vortexing.
3. The reconstituted solution can be stored at 2-8°C for up to 1 week.
4. For long-term storage, it is recommended to further dilute the solution with a carrier protein (such as 0.1% BSA, 10% FBS, or 5% HSA) to a concentration of no less than 10 µg/mL and aliquot for storage at -20°C to -80°C for 3 to 6 months. If serum-free experiments are required, a 5% trehalose solution can be used as a carrier instead. Avoid repeated freeze-thaw cycles.

Storage

Physical Appearance	Storage	Stability
Lyophilized powder	-20°C to -80°C	1 year
Reconstitution (initial)	2°C to 8°C	Less than 1 week
Reconstitution (after dilution)	-20°C to -80°C	3 to 6 months

Note

The reagent is only used in the field of scientific research, not suitable for clinical diagnosis or other purposes.