

## Sterile Water

Cat #: C-CYP399

Size: 1.5 mL

Shipping: Blue Ice

### Product Overview

This product is sterile ultrapure water, colorless and transparent, with no foreign substances.

The product can be used to resuspend lyophilized recombinant cytokines, and the recommended dissolution concentration is not less than 100 µg/mL and not more than 1 mg/mL. After adding an appropriate volume of sterile water to the recombinant cytokine product, it should be left at room temperature for at least 20 minutes to fully dissolve. Do not vortex or vigorously shake.

### Product Information

**Sterile:** Ready-to-use, sterilized by 0.22 µm filtration

**Application:** Reconstitution of lyophilized recombinant cytokines in suspension

**Storage:** Room temperature or 2-8°C, for 1 year

### Note

The reagent is only used in the field of scientific research, not suitable for clinical diagnosis or other purposes.